

RET[®] control-visc

The Magnetic Stirrer for Scientists

IKA[®]

شرکت ابزار طب

وارد کننده تجهیزات
آزمایشگاهی - تحلیلی - علمی
تاسیس ۱۳۶۳



تهران، خیابان مطهری، خیابان فجر، پلاک 38

021-88832100



021-88321727

021-88321728



www.abzarteb.com



[abzar.teb](https://www.instagram.com/abzar.teb)



Phil S. Baran
Ph.D.

Recipient of MacArthur
Genius Grant



Baran
uses IKA[®]
equipment



designed
to work perfectly


The RET® control-visc is the safest, strongest and most intelligent magnetic stirrer in its class.

The RET® control-visc is a magnetic stirrer whose remarkable technical functions have been developed for demanding applications. The unit mainly focuses on three core competences:
1. Safety, 2. Power, 3. Intelligence.

This is made possible by
> using high performance electronic components,
> intelligent heating technology,
> a motor designed specifically for a variety of applications (including high-viscous fluids) and
> high quality standards applied during the production process.

Insulated composite heating plate

With the unique structure of the composite heating plate, the RET® control-visc minimizes the loss through eddy currents when heating and stirring. The integrated high-tech insulation optimizes the heat transfer into the medium by minimizing thermal losses. The built-in heating foil ensures an even temperature allocation on the heating plate.



Sealed housing to protect motor and display




An integrated and patented weighing function allows the user to measure weight changes of up to 5,000 g



Torque trend measurement Viscosity changes in the medium can be measured by using a torque measurement device. The results can be depicted on the display



An RS 232 and USB interface enable connecting the unit to a PC for operating and updating the device



RET® control-visc with high-quality stainless-steel heating plate surface

The stainless steel surface of the composite plate enables the most efficient heat transfer to the medium and results in the fastest possible heating of the medium.



RET® control-visc white with ceramic coated heating plate

The RET® control-visc white offers a ceramic coated heating plate. The white surface helps to recognize color changes of fluids in a glass vessel.


شرکت ابزار طب
 گروه دانش تجهیزات
 پزشکی و دندانپزشکی
 تأسیس ۱۳۸۲
 تهران، خاندان مطهری، خاندان مهر، پلاک ۳۸
 ۰۲۱-۸۸۳۲۱۰۰
 ۰۲۱-۸۸۳۲۱۷۲
 ۰۲۱-۸۸۳۲۱۷۲
 www.abzar.teb.com
 abzar.teb

> The RET® control-visc offers excellent safety

The device comes with a coated and sealed housing which protects liquids from entering into the magnetic stirrer. Overheating is prevented by several integrated technical features. In the case of a malfunction, the device shuts down automatically and shows the error code on the TFT display. The integrated safety features also allow for an unsupervised operation of the RET® control-visc.



> Sealed housing



In case of a liquid overflow a built-in drainage protects the electronic components of the device.

- > Liquids cannot get inside the unit
- > Components are safe
- > Isolated drain




> Three temperature safety protection features

"Safety Temperature" is an adjustable temperature safety circuit that prevents from exceeding a specified set temperature. The safety temperature can be adjusted by using a special tool included in the product delivery.



"Set temperature" can be adjusted easily. It is used to safely heat the medium until the set temperature is reached.

"Overheating protection" Should the internal temperature of the RET® control-visc exceed the permissible temperature that would damage the internal electronic components, the heating power is reduced automatically.



> Safe operation

Operating modes

The unit is equipped with three operating modes:
A Mode: regular operation, all values can be directly changed
B Mode: all settings are stored when the device is switched off or loses power, functions are restored when the unit is powered ON again.
C Mode: if operating in C-Mode the set values are not changeable. When restarting the device these values are still fixed. In order to change the parameters, the software mode has to be changed to A or B through the display menu.



Password protection

Menu access can be password protected. If enabled, users cannot change any settings without password.

Adjustable limits

Limits can be set for speed and temperature. It is possible to set a minimum value for each parameter.

Lock button protects set parameters

YOUR BENEFITS

Coated and sealed housing

- > Liquids can not get inside the unit
- > Components are safe
- > Isolated drain
- > Protection class IP42

Three temperature safety protection features

- > Highest possible safety especially when working with easily flammable liquids
- > Manually adjustable safety circuit
- > Overheating protection for electronic components

C Mode advantages

- > Protected against changes to set values
- > Values are still fixed after restarting the device, suitable for serial testing
- > Automated restart after power outage to operating mode and set values



> RET® control-visc is the strongest magnetic stirrer in its class

Three components provide for an extraordinarily powerful magnetic stirrer:

1. high performance EC motor with 12W output
2. high performance internal transformer providing efficient power
3. composite heating plate with minimal eddy current losses

The unique structure of the insulated heating plate results in faster heating than other magnetic stirrers.

> Stirring performance
Powerful EC motor with high performance internal transformer



Compact and closed composite heating plate, combined with an advanced heating foil and engineered insulation, ensures an even temperature distribution on the heating plate.

- YOUR BENEFITS**
- Highly powerful and energy efficient
 - > High stirring speed stability
 - > Fast heating times
 - > High temperature stability
 - > Motor/transformator/composite heating plate = high performance of stirring and heating
 - > Engineered heating plate insulation
 - > Optimized heating through intelligent product design

> Heating performance
powerful and efficient heat transfer into the sample

Heating rate
 7 K/min for 1l H₂O at 600 W

شرکت ابزار طب
 راه توسعه تجهیزات
 آزمایشگاهی - تخصصی - صنعتی
 تأسیس ۱۳۶۳

تهران، خدایان مطهری، خدایان فجر، پلاک ۳۸
 021-88832100
 021-88321727
 021-88321728
 www.abzarreb.com
 abzar.teb

RET® control-visc | Intelligence



شرکت ابزار طب
 وزارت بهداشت، درمان و آموزش پزشکی
 تأسیس ۱۳۶۳
 تهران، خیابان مطهری، خدایان محور، پلاک ۳۸
 ۰۲۱-۸۸۸۳۲۱۰۰
 ۰۲۱-۸۸۳۲۱۷۲۷
 ۰۲۱-۸۸۳۲۱۷۲۸
 www.abzar.teb
 abzar.teb

> Easy operation with user-friendly display

The RET® control-visc continues the user-friendly tradition of operating the unit with two rotating knobs. They enable the easy and direct change of the most important parameters on the display menu.

The high-resolution display has easy to understand icons that allow for simple navigation through the menu, as well as allow for adjusting display settings, using the weighing or torque trend measurement functions, or changing the display language.

Easy operation and display of all relevant information at one glance



Set temperature

Set speed

Lock function Operating mode Temperature probe connected

Shows safety temperature and up to three user-defined parameters

Actual heating plate or probe temperature

Heating activated

Actual speed

Motor status

Safe Temp: 250 °C Timer: 00:03:58
 pH-Value: 7.0 Δ Torque: 100 %

233.8 actual °C 1700 actual rpm


240.0 set °C 1700 set rpm

> IKA® Magnetic Stirrer with Scientists

The Scripps Research Institute (TSRI) is a nonprofit research institution whose philosophy emphasizes the creation of basic knowledge in the biosciences for its application in medicine, the pursuit of fundamental scientific advances through inter-disciplinary programs and collaborators, and the education and training of researchers preparing to meet the scientific challenges of the future.

For more information on The Scripps Research Institute, see: www.scripps.edu

For more information on Phil S. Baran and the Baran Laboratory, see: www.scripps.edu/baran



The heart of the RET® control-visc is the ARM-based microcontroller which is also used in smart phones and tablets. The use of the ARM-based microcontroller technology provides the intelligence of simple navigation, firmware update possibility, weighing and torque trend measurement.



> Advanced technology

Integrated ARM-based microcontroller

The RET® control-visc uses technology which is used in smart phones or tablets. Two integrated ARM-based microcontroller along with a graphic controller are the base for all intelligent functions within the RET® control-visc. They provide for speed, energy efficiency and powerful performance. When selecting components, the IKA® engineers focus on quality, safety and reliability.

> Intelligent features



Integrated and patented weighing function

Performs simple weighing tasks without taking the sample off the device.

Measure weight changes of up to 5,000g	
tolerance	<500 g : +/-1g
	>500 g : +/-5g



Torque trend measurement

Relative viscosity changes can be measured with this feature by using a torque trend measurement device. Results can be depicted on the display. Useful for long term studies, test results can be documented through labworldsoft®. Reproducibility with max. deviation of +/-1%.



Stirring bar decoupling detection

The stirring function stops briefly when a decoupling occurs. It will automatically resume to the previously set speed when the stir bar is recoupled. Useful for long-term studies and when working with non-transparent fluids.



Firmware update tool

- > Keep your device up-to-date
- > Software upgrade features

The RET® control-visc has various interfaces: USB, RS 232 and Bluetooth



Firmware update tool



Interface

The RET® control-visc has an RS232 and USB interface; connect the unit to a PC for controlling and updating the device.



شرکت ابزار طب
آرام شده تجهیزات
برای آزمایشگاه، مطبها، مدارس
تاسیس ۱۳۶۳

تهران، خدایان مسپوری، خدایان فجر، پلاک 38

021-88832100
021-88321727
021-88321728
www.abzarteb.com
abzar.teb

YOUR BENEFITS

Intelligent solutions

- > User-friendly
- > Simple navigation and easy operation
- > Multilingual task menu
- > User-defined display settings
- > Integrated patented weighing function
- > Unique torque trend measurement
- > Stirring bar decoupling detection
- > Firmware update tool
- > labworldsoft® compatible



Integrated features



Safety



Power



Intelligence



- > **Sealed housing**
The housing is sealed to protect electronic components inside the device. Protection class IP 42.
- > **Three temperature safety protection features**
Set temperature for sample heating, adjustable safety circuit to avoid overheating of the heating plate, and overheating protection of internal components.
- > **Three different operating modes**
A Mode for regular operation, B Mode maintains preset parameters by the user from last operation when unit is powered on, C Mode for automatic restart after power outage to resume operation at preset settings.
- > **Password protection**
Menu access can be password protected. If enabled, users cannot change any settings without password.
- > **Lock function**
"Key" lock to prevent changes of set values.
- > **High-efficient motor power performance**
enables more stirring power (motor rating output: RH basic = 2 W, RCT basic = 9 W, RET* control-visc = 12 W)
- > **Innovative heating plate**
Insulated composite heating plate results in efficient heating of the sample with minimal eddy current losses.
- > **Heating plate with stainless steel surface**
enables quickest and safest heating of the sample.
- > **Temperature control with dual temperature sensor**
Simultaneous control of the temperature of the heating plate and the sample.
- > **High-resolution TFT display**
Lets the user see all relevant data clearly and simultaneously.
- > **Bluetooth interface**
Integrated Bluetooth interface (to use with Labworldsoft® and PC connection)
- > **Integrated weighing function**
Perform simple weighing tasks without taking the sample off the device.
- > **Torque trend measurement**
Resolve viscosity changes can be measured with this feature.
- > **Stirring bar decoupling detection**
The stirring function stops briefly when a decoupling occurs. It will automatically resume to the previously set speed when the stir bar is recoupled.
- > **pH measurement**
Measuring pH values is possible with pH sensors (BNC-connector); available through IKA® or other manufacturers.
- > **Labworldsoft®**
Specifically designed and developed by IKA®, this software allows for the RET* control-visc and other lab equipment from other manufacturers to be operated. For more information, please go to <http://www.ika.com>.
- > **Software update**
Keep your device software up-to-date with the integrated firmware update function. You can update through the USB interface of the stirrer and your PC.
- > **Multilingual menu**
The user can change the display language, there are 9 languages to choose from.
- > **Display view**
can be modified by the user, certain values can be shown or hidden.
- > **Error code display**
When an error occurs, the code is shown on the display. Please refer to the manual for further details.
- > **Reset**
The device can be reset to factory defaults.
- > **"HOT" warning**
When the device is off and the heating plate is hot, the display shows the warning "HOT" and the current heating plate temperature. The display turns off completely when the temperature of the heating plate drops below 50 °C.
- > **Interval mode**
The stirring function can be programmed to stop and automatically restart in intervals adjustable by the user.
- > **Sensor calibration**
Temperature and pH sensors can be calibrated through the RET* control-visc against a known calibrated source (i. e. temperature device and pH buffer solution respectively).
- > **Timer function**
Stop the heating process automatically after a specified time (max. 99:59:59)



شركت ابزار طب
توران، خندان مسوری، خندان تور، پلاک 38
021-88832100
021-88321727
021-88321728
www.abzarteb.com
abzar.teb

Technical data

Technical data	RET [®] control- <i>visc</i> / RET [®] control- <i>visc</i> white
Max. stirring quantity (R2)	20 l
Motor rating input / output	22 / 12 W
Speed range	50 – 1700 rpm
Heat output	600 W
Temperature range	RT – 340 °C
Adjustable safety circuit	50 – 380 °C
Control accuracy with sensor	PT 1000 ± 0.2 K
Heating plate material	stainless steel 1.4301 / white ceramic
Heating plate dimensions	Ø 175 mm
Dimensions (H x D x B)	160 x 270 x 85 mm
Weight	3 kg
Permissible ambient temperature	5 – 40 °C
Price	1.043,- EUR 1.066,- EUR

Protection class according to DIN EN 60529



3 Year warranty*



Protection class according to DIN EN 60529: IP 42

* 2+1 years after registering at www.ika.com/registration, gas/water and wiring parts included!

شرکت ابزار طب
تهران - خیابان مطهری - مقابل بیمارستان پارس
021-88832100
021-88321727
021-88321238
www.abzartob.com
Abzartob



- The RET[®] control-*visc* comes with the following accessories:
- PT 100.70 temperature sensor
Temperature sensor, stainless steel
 - Safety circuit tool
Tool to adjust safety circuit
 - H 104 protective cover (transparent)
Cover for RET[®] control-*visc*
 - Magnetic stirring bars
IKAION[®] 30 round PTFE coated
IKAION[®] 40 round PTFE coated
 - USB Cable
 - Power plug H 11
Mains cable, plug depending on region

Overview IKA[®] magnetic stirrers

IKA[®] offers a wide range of magnetic stirrers. Compare the following IKA[®] hotplate stirrers to help you find the most suitable unit for your application.



Technical data	RH basic 2	RH basic	RH digital
Display	Scale	Scale	LED
Max. heating plate temperature	320°C	320°C	320°C
Heat output	400 W	600 W	600 W
Max. stirring quantity (R2)	10 l	15 l	15 l
Heating plate material	stainless steel	stainless steel / composite / white ceramic	stainless steel / composite / white ceramic
Connection for ext. temp. sensor	–	ETS-D ✓	ETS-D x ✓
Integrated temperature regulation	–	–	–
Control accuracy with integrated temperature regulation	–	–	–
Interface for external control > bidirectional	–	–	–
Weighing, torque bend, pH	–	–	–
Firmware update tool	–	–	–
Protection class according to DIN EN 60529	IP 21	IP 21	IP 21

Basic stirring and heating functionality

Basic stirring and heating functionality with composite heating plate in stainless steel or white ceramic coated. External temperature control possible by connecting a contact thermometer (only ETS-D series).

Basic stirring and heating functionality with composite heating plate in stainless steel or white ceramic coated. External temperature control is possible by connecting a contact thermometer (only ETS-D series). Digital LED display for speed and temperature.



Technical data	C-MAG HS 7	C-MAG 7 digital	RCT basic	RET [®] basic	RET [®] control- <i>visc</i>
LED Scale	500°C	500°C	LED	LED	high res TFT
Max. heating plate temperature	320°C	320°C	320°C	340°C	340°C
Heat output	1000 W	1000 W	600 W	600 W	600 W
Max. stirring quantity (R2)	15 l	15 l	20 l	20 l	20 l
Heating plate material	white ceramic	white ceramic	aluminum alloy	stainless steel / composite	stainless steel / composite / white ceramic
ETS-D x ✓	PT 1000 ✓	PT 1000 ✓	PT 1000 ✓	PT 1000 ✓	PT 1000 ✓
precision	± 0.5 K ✓	± 0.5 K ✓	± 1 K ✓	± 1 K ✓	± 0.2 K ✓
ETS 2320 USB	–	–	–	–	–
Weight	–	–	–	–	–
Firmware update tool	–	–	–	–	–
Protection class according to DIN EN 60529	IP 21	IP 21	IP 42	IP 42	IP 42

Basic stirring and heating functionality with full-ceramic square plate to achieve higher temperatures. External temperature control is possible by connecting a contact thermometer with precise control accuracy (ETS-D series). Digital LED display for heating plate temperature.

Basic stirring and heating functionality with full-ceramic square plate to realize higher temperatures. External temperature control is possible by connecting the included temperature sensor (PT 1000). Digital LED display for speed and temperature.

Extended stirring and heating performance with aluminum alloy heating plate. External temperature control is possible by connecting the included temperature sensor (PT 1000). Digital LED display for speed and temperature.

Optimized stirring and extended heating performance with composite stainless steel heating plate. External temperature control is possible by connecting the included temperature sensor. Digital LED display for speed and temperature.

Optimized stirring and heating performance with composite heating plate in stainless steel or white ceramic coated. External temperature control is possible by connecting the included temperature sensor (PT 1000). The unit is capable of high-precision temperature control (± 0.2K). High-resolution TFT display enables easy operation.

Please visit www.ika.com for more information on IKA's magnetic stirrers and accessories.

Accessories



Sensors

	Ident. No.	Price
1 PT 1000.50*	0003367000	201,- EUR
Temperature sensor, dual, stainless steel, Ø 3 mm, 230 mm length		
2 PT 1000.51*	000337700	233,- EUR
Temperature sensor, dual, stainless steel, glass-coated, Ø 3 mm		
PT 1000.53*	000499800	243,- EUR
Temperature sensor, dual, coated with SafeCoat, 230 mm length		
PT 100.50	000201900	169,- EUR
Temperature sensor, stainless steel, d. 3 mm, 230 mm length		
3 PT 100.51*	000200300	201,- EUR
Temperature sensor, stainless steel, glass-coated, Ø 8 mm, 230 mm length		
PT 100.53	000499700	206,- EUR
Temperature sensor, stainless steel, coated with SafeCoat, Ø 3 mm, 230 mm length		
pH 5000.1	002000404	on request
pH sensor, pH 0-14, BNC connector, 1 m cable		



Stainless steel vessels

	Ident. No.	Price
1 H 1000	000444401	274,- EUR
Beaker, stainless steel, 1, inner Ø 160 mm, 79 mm height		
2 H 1500	000444501	265,- EUR
Beaker, stainless steel, 1.5 l, inner Ø 140 mm, 160 mm height		
3 H 3000	000444503	266,- EUR
Beaker, stainless steel, 3 l, inner Ø 180 mm, 119 mm height		
4 H 5000	000444505	329,- EUR
Beaker, stainless steel, 5 l, inner Ø 220 mm, 142 mm height		
5 H 8000	000444508	503,- EUR
Beaker, stainless steel, 8 l, inner Ø 265 mm, 162 mm height		

* Only available for RET* control-vice



Synthesis blocks

Dry heating block square series		Ident. No.	Price
H 135.10 Square carrier without handle	–	002500832	99.50 EUR
H 135.101 Block 16 x 4 ml	15.2 mm	002500626	37,- EUR
H 135.102 Block 16 x 8 ml	17.5 mm	002500627	40.50 EUR
H 135.103 Block 8 x 16 ml	20.5 mm	002500628	40.50 EUR
H 135.104 Block 4 x 20 ml	28.5 mm	002500629	34.50 EUR
H 135.105 Block 4 x 20 ml	28.5 mm	002500630	40.50 EUR
H 135.106 Block 4 x 40 ml	28.5 mm	002500631	45,- EUR
H 135.107 Block 130 ml	–	002500632	36,- EUR
H 135.108 Block 250 ml	–	002500633	42.50 EUR

100 ml Flask heating block series

	Max. Outer Ø	Inner Ø	Ident. No.	Price
H 135.20 Flask carrier 100 ml without handle	142 mm	–	002500634	87,- EUR
H 135.201 Flask inlet 10 ml	–	33.8 mm	002500636	20.50 EUR
H 135.202 Flask inlet 25 ml	–	43.8 mm	002500637	20.50 EUR
H 135.203 Flask inlet 50 ml	–	52.8 mm	002500638	20.50 EUR

500 ml Flask heating block series

H 135.30 Flask carrier 500 ml without handle	142 mm	–	002500639	99.50 EUR
H 135.301 Flask inlet 100 ml	–	66.3 mm	002500641	46.50 EUR
H 135.302 Flask inlet 250 ml	–	88 mm	002500642	47.50 EUR

1000 ml Flask heating block series

H 135.40 Flask carrier 1000 ml without handle	166.3 mm	–	002500633	144,- EUR
H 135.401 Flask inlet 500 ml	–	108 mm	002500644	81.50 EUR

2000 ml Flask heating block series

H 135.50 Flask carrier 2000 ml without handle	194.7 mm	–	002500634	178,- EUR
H 135.501 Flask inlet 1000 ml	–	131.8 mm	002500645	119,- EUR

3000 ml Flask heating block series

H 135.60 Flask carrier 3000 ml without handle	218.3 mm	–	002500635	256,- EUR
H 135.601 Flask inlet 2000 ml	–	162 mm	002500648	189,- EUR

5000 ml Flask heating block series

H 135.70 Flask carrier 5000 ml without handle	299.7 mm	–	002500636	389,- EUR
H 135.701 Flask inlet 3000 ml	–	194.5 mm	002500651	331,- EUR



Slitting bars

IKAFLOW® round	Ident. No.	Dimensions (L x Ø)	Price
IKAFLOW® 10 Set (5 Pcs) round PTFE-coated	000448800	10 x 6 mm	6.40 EUR
IKAFLOW® 15 Set (5 Pcs) round PTFE-coated	000448870	15 x 6 mm	6.40 EUR
IKAFLOW® 20 Set (5 Pcs) round PTFE-coated	000448900	20 x 8 mm	6.40 EUR
IKAFLOW® 25 Set (5 Pcs) round PTFE-coated	000448960	25 x 8 mm	9.70 EUR
IKAFLOW® 30 Set (5 Pcs) round PTFE-coated	000449000	30 x 8 mm	8.10 EUR
IKAFLOW® 40 Set (5 Pcs) round PTFE-coated	000449100	40 x 8 mm	9.80 EUR
IKAFLOW® 50 Set (5 Pcs) round PTFE-coated	000449200	50 x 8 mm	11.70 EUR
IKAFLOW® 80 Set (5 Pcs) round PTFE-coated	000449300	80 x 10 mm	26,- EUR
IKAFLOW® glass 25 Set (5 Pcs) round	000449220	25 x 8 mm	11.40 EUR
IKAFLOW® glass 30 Set (5 Pcs) round	000449240	30 x 8 mm	11.90 EUR
IKAFLOW® glass 40 Set (5 Pcs) round	000449300	40 x 8 mm	13,- EUR
IKAFLOW® glass 50 Set (5 Pcs) round	000449300	50 x 8 mm	14,- EUR
IKAFLOW® 15 Set (5 Pcs) power Sema PTFE-coated	000449100	15 x 9 mm	18.60 EUR
IKAFLOW® 30 Set (5 Pcs) power Sema PTFE-coated	000449320	30 x 12 mm	28,- EUR
IKAFLOW® 80 Set (5 Pcs) power Sema PTFE-coated	000449400	50 x 21 mm	34,- EUR

IKAFLOW® Slide round

IKAFLOW® 10 Set (5 Pcs) slide round PTFE-coated	000449300	25 x 6 mm	8,- EUR
IKAFLOW® 15 Set (5 Pcs) slide round PTFE-coated	000449400	30 x 6 mm	8.60 EUR
IKAFLOW® 20 Set (5 Pcs) slide round PTFE-coated	000449420	40 x 8 mm	14.40 EUR
IKAFLOW® 30 Set (5 Pcs) slide round PTFE-coated	000449440	50 x 8 mm	16.60 EUR

IKAFLOW® ellipse

IKAFLOW® 10 Set (5 Pcs) ellipse PTFE-coated	000449400	20 x 10 mm	12.70 EUR
IKAFLOW® 15 Set (5 Pcs) ellipse PTFE-coated	000449500	25 x 12 mm	13.90 EUR
IKAFLOW® 20 Set (5 Pcs) ellipse PTFE-coated	000449600	30 x 15 mm	18.30 EUR
IKAFLOW® 25 Set (5 Pcs) ellipse PTFE-coated	000449700	35 x 15 mm	23.50 EUR
IKAFLOW® 40 Set (5 Pcs) ellipse PTFE-coated	000449800	40 x 20 mm	29.50 EUR
IKAFLOW® 50 Set (5 Pcs) ellipse PTFE-coated	000449900	50 x 20 mm	41.50 EUR
IKAFLOW® 70 Set (5 Pcs) ellipse PTFE-coated	000449980	64 x 20 mm	60.50 EUR
IKAFLOW® 80 Set (5 Pcs) ellipse PTFE-coated	000449600	70 x 20 mm	73,- EUR

IKAFLOW® cross

IKAFLOW® 10 Set (5 Pcs) cross PTFE-coated	000449620	10 x 10 mm	10.90 EUR
IKAFLOW® 15 Set (5 Pcs) cross PTFE-coated	000449640	20 x 20 mm	13.90 EUR
IKAFLOW® 20 Set (5 Pcs) cross PTFE-coated	000449660	25 x 25 mm	15.90 EUR
IKAFLOW® 30 Set (5 Pcs) cross PTFE-coated	000449680	30 x 30 mm	20.50 EUR
IKAFLOW® 38 Set (5 Pcs) cross PTFE-coated	000449700	38 x 38 mm	25.50 EUR

IKAFLOW® bone

IKAFLOW® 17 Set (5 Pcs) bone PTFE-coated	000449720	37 x 8 x 20 mm	37.50 EUR
IKAFLOW® 54 Set (5 Pcs) bone PTFE-coated	000449740	54 x 8 x 20 mm	45,- EUR

TRIKA®

TRIKA® 25 Set (5 Pcs) PTFE-coated	000449700	25 x 8 mm	6.70 EUR
TRIKA® 40 Set (5 Pcs) PTFE-coated	000449800	35 x 9 mm	13.50 EUR
TRIKA® 55 Set (5 Pcs) PTFE-coated	000449900	55 x 14 mm	21.50 EUR
TRIKA® 80 Set (5 Pcs) PTFE-coated	000449960	80 x 18 mm	30,- EUR

IKAFLOW® beaker

	Ident. No.	Dimensions (Ø x H)	Price
IKAFLOW® 67 beaker PTFE-coated	000449700	67 x 21 mm	208,- EUR
IKAFLOW® 74 beaker PTFE-coated	000449700	74 x 29 mm	215,- EUR
IKAFLOW® 103 beaker PTFE-coated	000449800	103 x 32 mm	310,- EUR
IKAFLOW® 125 beaker PTFE-coated	000449820	125 x 48 mm	396,- EUR

IKAFLOW® crown

IKAFLOW® 9 Set (5 Pcs) crown PTFE-coated	000449840	9 x 6 mm	9,- EUR
--	-----------	----------	----------------



Lifts

LABLIFT m	0004022400	158,- EUR
Manual scoop lift		



Other accessories

	Ident. No.	Price
H 16 V	0001545100	47.50 EUR
Support rod, Ø 10 mm, 450 mm length		
H 16.1	0005005000	12.80 EUR
Extension for support rod		
H 44	0002437700	21,- EUR
Boss head clamp		
H 38	0003547700	21,- EUR
Holding rod		

YOUR BENEFITS

- Dual sensors (PT 1000.50, PT 1000.51, PT 1000.53) can be connected
 - > Precise temperature control
 - > Safe handling of highly temperature sensitive fluids
- IKA stainless steel vessels ensure:
 - > Optimized heat transfer
 - > Outstanding magnetic adhesion and
 - > Fits securely on the heating plate




IKA® offers more

شرکت ابزار طب
 وارد کننده تجهیزات
 آزمایشگاهی، تجهیزات - علمی
 تأسیس ۱۳۶۳

تهران، خیابان مطهری، خیابان فجر، پلاک 38

021-88832100
 021-88321727
 021-88321728

www.abzarteb.com
 abzar.teb



Modern manufacturing

During manufacturing, IKA® focuses on high quality, not only with well-trained and experienced personnel, but also with standardized processes and quality checks.

The assembly of the printed circuit boards is fully automated and includes an automated 100% quality control check of every PCB.

Customizing Center

It is important that IKA products work for your application. We are introducing a new program: product solutions tailored to your needs.

Should you not find the appropriate device in our standard product range, please send us your requested specifications through the online form. Our team will determine its feasibility and offer a solution to you.

Please visit www.ika.com/customizingcenter to review already implemented product modifications.

Worldwide service network – direct contact in your region

Our dedicated team of engineers provides comprehensive worldwide technical service. Please feel free to contact IKA® directly or your dealer in case of any service questions.

For spare parts IKA® guarantees 10 years of availability. In the event of an equipment malfunction or technical questions regarding devices, maintenance and spare parts, please call us at 00 8000 4524357 (00 8000 IKAHELP) or send an email to service@ika.com

IKA® Application Support

Our Application Center spans 400sqm and offers modern facilities for presenting and testing lab devices and processes. This brings us even closer to our customers and improves our service. Here, prospective buyers and customers can test processes that involve stirring, shaking, dispersing, grinding, heating, analyzing and distilling.

Call us at 00 8000 4522777 (00 8000 IKAAPPS) or send an email to applicationsupport@ika.com or visit our website at www.ika.com/applicationsupport

Prices valid until December 31st, 2014
Subject to technical changes and
alteration of prices

IKA®+

Ordering made easy!
For more information about
our products and to place
your order, please visit:
www.ika.com

00310



شرکت ابزار طب
وارد کننده تجهیزات
آزمایشگاهی - تحقیقاتی - علمی
تاسیس ۱۳۶۳

تهران، خیابان مطهری، خیابان فجر، پلاک 38



021-88832100

021-88321727

021-88321728



www.abzarteb.com



[abzar.teb](https://www.instagram.com/abzar.teb)



www.facebook.com/IKAworldwide

IKA®-Werke GmbH & Co. KG
Janke & Kunkel-Str. 10
79219 Staufen
Germany

Tel. +49 7633 831-0
Fax +49 7633 831-98

sales@ika.de
www.ika.com

IKA® German technology

