

## NSZ-606 Zoom Stereo Microscope



شرکت ابزار طب  
وارد کننده تجهیزات  
آزمایشگاهی - تحقیقاتی - علمی  
تاسیس ۱۳۶۳



تهران، خیابان مطهری، خیابان فجر، پلاک 38



021-88832100

021-88321727

021-88321728



[www.abzarteb.com](http://www.abzarteb.com)



[abzar.teb](https://www.instagram.com/abzar.teb)

## Description

- With high quality optical components, providing sharp and striking contrast image and ensure flat image at big depth of field, ease the tiredness with ergonomic design.
- With LED light for both incident and transmitted illumination, Providing even illumination and life expectancy can reach 6000 hours.
- With Complete eyepieces, objectives and accessories, it is an ideal instrument for industry assembling, inspection and teaching field.

## Specification:

		NSZ-606
<b>Drawing Tube</b>	Binocular Viewing Head, Inclined at 45°	•
	Trinocular Viewing Head, Inclined at 45°	○
<b>Eyepiece</b>	Extra wide field eyepiece EW10×/22	•
	WF15×	○
	WF20×	○
<b>Zoom Objective</b>	0.8×-5×	•
<b>Zoom Ratio</b>	6.3 : 1	•
<b>Working Distance</b>	115mm	•
<b>Interpupillary Distance</b>	57-75mm	•
<b>Auxiliary objective</b>	0.5×	○
	2×	○
<b>Glass Insert</b>	Glass Insert, Diameter 100mm	•
<b>Illumination</b>	Incident Illumination 100V-240V/ LED	•
	Transmitted Illumination 100V-240V/ LED	•
<b>Photo Attachment</b>		○
<b>Video Attachment</b>		○

Note: • Standard Outfit, ○ Optional



تهران، خیابان مطهری، خیابان فجر، پلاک 38



021-88832100

021-88321727

021-88321728



www.abzarteb.com



abzar.teb